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Target (Analyte)	Sensitivity	Specificity
Bacteria & Bacterial Toxins	97.7%	99.9%
Viruses	97.7%	99.9%
Campylobacter	97.7%	99.9%
Escherichia coli O157	100%	100%
Shiga-like Toxin	99.2%	99.9%

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
Target (Analyte)	Sensitivity	Specificity
Bacteria & Bacterial Toxins		
Shigella	97.7%	97.8%
Campylobacter	97.5%	97.8%
Verotoxin-producing E. coli (VTEC)	N/A*	100.0%
Enterohaemorrhagic E. coli (EHEC) (EHEC)	N/A*	97.2%
E. coli O157	99.2%	99.8%
Shiga-like Toxin-producing E. coli (STEC) (sh1/hly2)	100.0%	99%
Salmonella	92.1%	98.1%
Clostridium difficile/toxin A/B	97.7%	94.9%
Vibrio cholerae	N/A*	100.0%
Toxins		
Adenovirus 40/41	100.0%	100.0%
Rotavirus A	94.7%	99.8%
Norovirus GII/IV	93.5%	99%
Parasites		
Giardia lamblia	100.0%	99.2%
Cyclospora	97.5%	100%
Enterococcus faecalis	100.0%	99.4%

*Due to the sample size, clinical sensitivity was not assessed for ETEC, Verotoxin-producing E. coli and Vibrio cholerae. However, analytical accuracy for these analytes was demonstrated in the limit of detection and specificity studies with cultural isolation or plasmids.

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

- [1] World Health Organization. Diarrhoeal Disease, Fact sheet N°330. August 2009. Available at: <http://www.who.int/media-centre/factsheets/fs330/en/index.html>
- [2] U.S. Department of Health and Human Services. 2010 National Statistics on Intestinal Infections. Available from HCUPnet at: <http://hcupnet.ahrq.gov/HCUPnet.jsp>
- [3] Eric C. Claas, Carey-Ann D. Burnham, Tony Mazzulli, Kate Templeton, and Francois Topin, 2013. Performance of the xTAG® gastrointestinal pathogen panel, a multiplex molecular assay for simultaneous detection of bacterial, viral, and parasitic causes of infectious gastroenteritis. *J. Microbiol. Biotechnol.* (2013), 23(7); 1041–1045
- [4] J.M. Bible, J. S. Ramble, L. Carlton, Carew, P. Cliff, E.P. Halligan. An evaluation of the Luminex xTAG® Gastrointestinal Pathogen Panel at a London teaching hospital 2011: The comparative performance of a rapid molecular multiplex assay and current standard laboratory investigations of gastroenteritis. Available at: <http://www.luminexcorp.com/Products/Assays/ClinicalDiagnostics/xTAGGPP/>
- [5] Mengelle C, 55Mansuy JM, 55Prere MF, 55Grouteau E, 55Claudet I, 55Kamar N, 55Huynh A, 55Plat G, 55Benard M, 55Marty N, 55Valentin A, 55Berry A, 55Izopet J. Simultaneous detection of gastrointestinal pathogens with a multiplex Luminex-based molecular assay in stool samples from diarrhoeic patients. *Clin Microbiol Infect.* 2013 Apr 29.

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